GLUCOSE ISOMERASE from STREPTOMYCES OLIVACEUS

Prepared at the 28th JECFA (1984), published in FNP 31/2 (1984) and in FNP 52 (1992). An ADI 'acceptable' was established at the 29th JECFA (1985).

SYNONYMS	Xylose isomerase
SOURCES	Produced by the controlled fermentation of Streptomyces olivaceus
Active principles	Xylose isomerase (glucose isomerase)
Systematic names and numbers	D-Xylose ketol-isomerase (EC 5.3.1.5)
Reactions catalyzed	D-Xylose and D-glucose are converted to D-xylulose and D-fructose, respectively
DESCRIPTION	Off-white to brown granules (immobilized preparation) or liquids, insoluble in water (granules), ethanol, chloroform and ether The immobilized preparations are rendered insoluble in water by treatment with gelatine (carrier) and glutaraldehyde (immobilization agent).
FUNCTIONAL USES	Enzyme preparation Used in the preparation of high fructose corn syrup and other fructose starch syrups.
GENERAL SPECIFICATIONS	Must conform to the <i>General Specifications for Enzyme Preparations used in Food Processing</i> (see Volume Introduction)
CHARACTERISTICS	
IDENTIFICATION	
Glucose isomerase activity (Vol. 4)	The sample shows glucose isomerase activity
PURITY	
<u>Glutaraldehyde</u> (Vol. 4)	Passes Limit Test for Glutaraldehyde from Immobilized Glucose Isomerases crosslinked with Glutaraldehyde